

Checklist

Delivery conditions for reclining ergometers

Dear customer,

To help ensure straightforward delivery of your new reclining ergometer, we need some information from you.

Please complete and return this checklist with accurate information. Missing or incorrect information may lead to unforeseen problems during the delivery or to exceeding the agreed delivery time, which can result in extra charges.

Please contact our sales team with any questions you may have regarding the delivery.

Tel.: +49-(0)-7431 - 9894 - 162

Fax: +49-(0)-7431 - 9894 - 128

E-mail: vertrieb@ergoline.com

Equipment will not be released for shipping until we receive your completed checklist!

1 Customer details

Name:	
Facility:	
Delivery address:	
Building:	
Postal code/town:	
Contact person:	
Telephone:	
E-mail:	
Order number:	

2 Access – unloading point

	Yes	No
Is the access road suitable for 20-ton trucks (12 m long, 3 m wide, 4 m high)?	<input type="checkbox"/>	<input type="checkbox"/>
Space for a trailer (12 m long)?	<input type="checkbox"/>	<input type="checkbox"/>
Unloading point located in a pedestrian zone?	<input type="checkbox"/>	<input type="checkbox"/>
Unloading point located in a construction site?	<input type="checkbox"/>	<input type="checkbox"/>
Unloading point located in a residential area?	<input type="checkbox"/>	<input type="checkbox"/>
Parking space available for trucks? (cannot be a no-parking zone!)	<input type="checkbox"/>	<input type="checkbox"/>
Passage through gateways/roofed structures?	<input type="checkbox"/>	<input type="checkbox"/>
Height of gateway/roofed structure?		m
Forklift/pallet truck available?	<input type="checkbox"/>	<input type="checkbox"/>

Delivery conditions for reclining ergometers

Characteristics of unloading point?		
Characteristics of transport paths?		
Do you have any on-site images of the access, unloading point, transport paths and set-up location? (If yes, please send the images together with the completed checklist to ergoline GmbH)	Yes <input type="checkbox"/>	No <input type="checkbox"/>

3 Local conditions (outside)

Ground floor	Yes <input type="checkbox"/>	No <input type="checkbox"/>			
Landing	Yes <input type="checkbox"/>	No <input type="checkbox"/>	Height:	cm	
Staircase	Yes <input type="checkbox"/>	No <input type="checkbox"/>	Number of steps:		
Type of staircase (see Item 7)			Stairs lead	Up <input type="checkbox"/>	Down <input type="checkbox"/>
Stair railing?	Yes <input type="checkbox"/>	No <input type="checkbox"/>			
Can the truck's lifting platform reach beyond the stairs?	Yes <input type="checkbox"/>	No <input type="checkbox"/>	Length of transport path:	m	
Covering required for outside transport path?	<input type="checkbox"/> Yes, meters		No <input type="checkbox"/>		
Does the floor have a sufficient load rating?	Yes <input type="checkbox"/>	No <input type="checkbox"/>			
Delivery opening	Height:	m	Smallest clearance height:	m	
	Width:	m	Smallest clearance width:	m	
Ramp	Yes <input type="checkbox"/>	No <input type="checkbox"/>	Raised platform	Yes <input type="checkbox"/>	No <input type="checkbox"/>

4 Local conditions (inside)

Ground floor	Yes <input type="checkbox"/>	No <input type="checkbox"/>	Floor level no.:		
Landing	Yes <input type="checkbox"/>	No <input type="checkbox"/>	Height:	cm	
Staircase	Yes <input type="checkbox"/>	No <input type="checkbox"/>	Number of steps:		
Type of staircase (see Item 7)			Stairs lead	Up <input type="checkbox"/>	Down <input type="checkbox"/>
Stair railing?	Yes <input type="checkbox"/>	No <input type="checkbox"/>	Raised platform	Yes <input type="checkbox"/>	No <input type="checkbox"/>
Length of transport path:	m				
Does the floor have a sufficient load rating?	Yes <input type="checkbox"/>	No <input type="checkbox"/>			
Covering required for outside transport path?	<input type="checkbox"/> Yes, meters		No <input type="checkbox"/>		
Delivery opening	Height:	m	Smallest clearance height:	m	
	Width:	m	Smallest clearance width:	m	

5 Elevator









Does the building have an elevator?	Yes <input type="checkbox"/>	No <input type="checkbox"/>			
Load capacity?	kg				
Height/width of elevator door?	Height:	m	Width:	m	
Height/width/depth of elevator cabin?	Height:	m	Width:	m	Depth: m

Delivery conditions for reclining ergometers

6 Set-up location

Is the room prepared?	Yes <input type="checkbox"/>	No <input type="checkbox"/>		
Is there turning space in front of the entrance to the set-up location? (min. 2.5 m x 2.5 m)	Yes <input type="checkbox"/>	No <input type="checkbox"/>	Yes <input type="checkbox"/>	No <input type="checkbox"/>
Building / floor / room				
Ground floor	Yes <input type="checkbox"/>	No <input type="checkbox"/>		
Landing	Yes <input type="checkbox"/>	No <input type="checkbox"/>	Height: cm	
Length of transport path:	m			
Does the floor have a sufficient load rating?	Yes <input type="checkbox"/>	No <input type="checkbox"/>	Floor covering:	
Covering required for transport path?	<input type="checkbox"/> Yes, meters No <input type="checkbox"/>			
Delivery opening	Height: m	Smallest clearance height:		m
	Width: m	Smallest clearance width:		m

7 Type of staircase (see Items 3 and 4) Important!

 <p>Type 1</p>	 <p>Type 2</p>
 <p>Type 3</p>	 <p>Type 4</p>
 <p>Type 5</p>	 <p>Type 6</p>
 <p>Type 7</p>	 <p>Type 8</p>

Date Signature/stamp