

### Patient Height/Saddle Types ergoselect 10/12

The standard scope of delivery of the ergoselect 10/12 includes the combination "Narrow saddle – standard plug-in holder" The following saddle combinations can be ordered as accessories:

Wide saddle – standard plug-in holder (erg705.979)

Wide saddle – extra short plug-in holder (erg705.985)

Narrow saddle – extra short plug-in holder (erg705.883)

All saddle types can be inserted in the saddle guide rail by means of the plug-in holder.

Due to the shape of the plug-in holder and horizontal displacement of the saddle, it is always possible to set the optimal supine position for the patient with respect to the pedal unit.

The angle between the axis formed by the upper body and the thigh should be about 10° when the pedal is in the bottom position.

The following leg lengths were determined by measurements (distance a in the illustration):

Wide saddle – standard plug-in holder: a = 77 cm

Wide saddle – extra short plug-in holder: a = 64 cm

Narrow saddle – extra short plug-in holder: a = 64 cm

### Important Information



The decisive factor for the patient's optimal supine position is the leg length (not the body height!). In combination with the optimal saddle, a wide range of body heights can be accommodated.

